



Name	PEARL SHELL SEALING MACHINE
Model	YSSM01A
Power Supply	220 / 240 V, 50 / 60 Hz
Power	500 W
Max. width of sealing	300 mm
Sealing band height	12 mm
Temperature range	90°C to 210°C
Product size	517x283x205 mm
Net weight	8.8 KGs



We focus on sterilization

NINGBO HAISHU YESON MEDICAL DEVICE CO., LTD
Add: No. 1 Huanzhen Dadao, Hengjie, Haishu 315181 Ningbo, China
Web: www.yeson-medicine.com Email: sale01@yeson-medicine.com
Cell/Whatsapp/Wechat: +86 13566025128




Provided by 



PEARL SHELL


sealing machine

Precisely made for dental clinics




Wide sealing width

Even if it's a 300 mm wide sterile pouch, it's easy to handle.



Sterile pouch free from loading

Thoughtful loading design of pouch rolls, free nurses from tedious and laborious work.



Groundbreaking paper-feeding

Eliminate the troubles during feeding paper, which also reduce the working load of clinical nurses. Efficiency are achieved at least 2-times faster.



Friendly human computer interface

LCD display, acoustic & visual prompt, paper-cutting counts, over-temperature protection, Standby protection, Overtime reminder.....



Joint Development of a **China-Austria** Team

